




**AUTOMATED SOLUTION FOR LIQUID  
HANDLING & CENTRIFUGAL PARTITION  
CHROMATOGRAPHY**

GILSON CPC

# Centrifugal Partition Chromatography (CPC) or Counter-Current Chromatography (CCC)

A CPC system is a chromatographic column with inlet and outlet as HPLC or flash chromatography. Main difference with classical LC: The stationary phase is liquid and non miscible to the eluent. To maintain the liquid stationary phase a CPC column must generate a centrifugal field and get cells with special design. A CPC column could be connect to prep HPLC systems.



PipetmaX

## Build Your Lab Assistant to Suit Your Needs.

A variety of tube racks, plates, blocks, pipetting heads, tips, and tip accessories make PIPETMAX® completely customizable to your application.

- PIPETMAX Pipette Heads, Tube Racks, Freezer Blocks, Tip Accessories, On-Bed Devices (Circulating Temperature Block and Adapters), Orbital Shaker, etc.

PIPETMAX is a pipette operator just like you. It can perform a variety of applications for you while you focus on more important things.

- qPCR
- Cell assays
- Enzymatic assays
- NGS library prep
- General liquid handling

